

# Assembly Instructions






**WARNING:** Please read all instructions before adjusting or using the walker, failure may result in injury or damage.  
 Maximum loading weight of seat: 100kg  
 Maximum loading weight of bag: 5kg

## Safety Instructions

1. Please make sure the walker is fully opened, seat is down and the folding hinge under seat is locked in the down position.
2. Be sure that both knobs are fully tightened and at proper height. Test walker for stability before applying full weight to hand grips.
3. Test the brake locks before operating the walker. When wheels are locked, brake shoe on back wheels should prevent wheels from turning (see operating instructions for brakes below). When brakes are released, the wheels should turn freely.
4. Always engage brake locks before sitting and do not use seat when unit is on slope or uneven ground.
5. Do not use the walker to walk backwards, down gradients or to climb stairs, curbs, or to go over obstacles. Serious risk of injury could occur.
6. The backrest is intended to provide back support while seated. The backrest is not intended to support full body weight.
7. All wheels must have contact with the floor and/or ground at all times.
8. Excess weight capacity can cause serious injury.

**NOTE:** Do not have anyone push you while you are seated on the walker. This is a walking aid only and is not to be used as a transportation device.

## Hardware List

A		1pcs	C		2pcs
			D		2pcs
B		1pcs	E		2pcs

## Part List



# Assembly Instructions



**Step 1:** Remove all components from carton: backrest, handle with tube, walker frame and basket.



**Step 2:** Open frame(Part A) and put down the folding hinge. Insert handgrips into the frame. Adjust the height of the handlebars for proper fit. To determine proper height, stand upright behind the rollator with shoulders relaxed, arms hanging at your side and eyes looking straight ahead. Standing in this position, the handgrips of the walker should be adjusted to a height even with the wrists. This will position the arms at an approximate 20 degree to 30 degree bend when using the walker.



**Step 3:** Align one of the holes to the hole in the frame and insert the bolt(D) to hold the handgrips assembly in place.



**Step 4:** Slide the washer(E) onto the other end of the bolt with the curved side facing the frame.



**Step 5:** Screw on the five point knob (C).



**Step 6:** Insert the backrest (B) so that the curved portion faces the front of the unit. Make sure that the button snaps firmly into the hole so the back is secure.



**Step 7:**

### Operating Instructions for brakes



1) To slow down, squeeze the brake handles up.


# Operating Instructions for brakes



2) To lock the brakes, push down on the brake handles until you feel or hear the brakes click into the locked position.



3) To release the brakes from the locked position, simply pull up on the brake handles.



4) Hold the thumbscrew firmly in place, do not turn.  
To tighten the brake, turn the adjustment nut clockwise;  
To loosen the brake, turn the adjustment nut counter-clockwise.

## Using the walker as a seat:

- 1) Push the walker to the desired location.
- 2) **LOCK** the brakes **BEFORE** sitting.

## Using the walker as a walker:

For proper positioning and use of the walker, imagine a line drawn across the back of the rear wheels. Part of one foot should be slightly ahead of that line as you walk.