front back

anko

Build & Paint Your Own Dinosaurs



TIPS:

- Don't disassemble all the wooden pieces at once or it will be difficult to match the pieces together for each dinosaur.
- 2. Use the included glue to hold the wooden pieces together when required.
- 3. Add some water to the paints if they look too dry or sticky.
- 4. Let the painted pieces dry completely after you have finished painting.
- 5. Keep the box nearby so that you can look at the picture for inspiration.
- 6. Freely mix up colour to create new ones.
- 7. It is recommended to paint lighter colours first and darker colours last.

DIRECTIONS 1:

- 1. Find a flat surface that you can use as your work area.
- 2. Cover work area with newspapers or paper towels,
- 3. Lay out all of your material in front of you.
- 4. Put on an apron before you paint.
- 5. Keep a cup of water to clean the paint brush after applying each colour.
- 6. Close the paint lid after use.


DIRECTIONS 2:

How to paint the dinosaurs:

paint the pieces on the board with the brush, you can refer to the picture on the packaging or create your own design.

2 How to build the dinosaurs:

Assemble the wooden pieces by matching the same numbers to each other, then repeat this action until the dinosaur is assembled.





148*260mm

anko

Build & Paint Your Own Dinosaurs







Note: Mixing Red and Yellow will make a close shade to Orange.

Mixing Blue, White and Yellow will make a close shade to Mint.

MARNING: CHOKING HAZARD SMALL PARTS, NOT FOR CHILDREN UNDER 3 YEARS.

ATTENTION: ALL ART MATERIALS LIKE PAINT, MAY CAUSE STAINS, ALMAYS PROTECT CLOTHES, CASPT, MORK SURVEY, ENHANCES CHOICE, CASPT, MORK SURVEY, ENHANCES CONTROL OF REVENT STAINS, USE A PROTECTIVE COVER TO PREVENT STAINS, PRODUCT MAY USE TUDINTLY PROM MADE SHOWN, PLEASE KIETS THE PACKAGINES FOR PUTUPER RETREMENT.

MADE IN CHINA MIND IS A TRADEMURIK OF MARTI AUSTRALIA UMITED.

CAUTION PRINCE PITS INMEDIATELY IF THE PAINT COMES INTO CONTICT WITH THEM, IF BREVIRTON PERSISTS, SEEN MEDICAL ATTENTION.

148*260mm