


8 x 9.5cm

# Lockable Journal

3+  
years


## OPEN JOURNAL DIRECTIONS:

1. Press the pass code into the front buttons at front of the journal (pass code is listed on the back of the journal).
2. Press the "confirm" button to release the lock, gently pull on the back cover to open.



## CLOSE JOURNAL DIRECTIONS:

1. Keep the pass code pressed, gently close the cover until the pin clicks into the keyhole.
2. Don't forget to reset all buttons to cover your pass code.



PLEASE NOTE THE PASS CODE CAN NOT BE CHANGED